



Patient Warmer

PW-10



Main Features

- Disposable blanket to keep off cross-infection
- Inflatable medical warming blanket can avoid partial intense heat
- Soft and collapsible, with an even array of air vents
- Two ways of usage, under and on the body
- Two size, full-length and half-length
- Hygienic blanket as lasting bacteria-proof filter inside warming equipment
- Three wind speed for adjustment
- The adoption of PTC heating technique, fast and effective
- Double self-protection from error and overheating, safe and sound
- Body Temperature can be detected and reflected on screen in real time
- Accumulated warming time display
- Three hotkeys of wind temperature for easy operation, 38°C, 40°C, 43°C, with the increment of

Technical Parameters

AC power supply	220V±10%.
50Hz±1Hz	
Power	1600AV
Environment temperature	5°C~40°C
Wind temperature adjustment range	20°C~50°C
Relative humidity	<80%
Atmospheric pressure	86KPa~106KPa
Wind volume	600dm ³ */min
Stretched length	2m



ACER BIOMEDICALS LLP

Office: Industrial Unit-80, Block-G, Sector-63, Noida, UP - 201301 (India)
Ph.: +91-120-4166921, E-mail: info@acerbiomed.com
Web.: www.acerbiomedicals.com, www.acerbiomed.com

ServiceHelpline
18008894232

