



Automatic Chemistry Analyzer

AC-120



Product features

Mini Fully Auto Chemistry Analyzer

Discrete/random access

STAT sample priority

Easy to use

User friendly interface

No computer

Internal thermal printer

Direct print report

All-in-one machine

Automatic Washing


No Washing (Option : SL120 D)

Internal printer

One button to test

Touch screen

Internal operating system



Working Conditions

Temperature: 10°C ~ 30°C

Humidity: ≤ 85%

Power Supply: 100~240V 50/60HZ

Water Consumption: 5L

Reaction/Reagent/Sample Handling :

- 81 cuvettes
- 24+5(Tube 12*100) sample positions
- 1 calibration position / 1 QC position / 1 ST position
- 54 positions for reagent
- 27 bottles with volume of 8ml
- 27 bottles with volume of 24ml
- LED light source with 8 channel transmissometry, Rear Spectrophotometry
- Extra 1 channel scattering (Option : SL120 C)
- Liquid level detection / Collision protection
- High pressure water flush interior the probe

Other

- Test speed : Constant 120T/H
- Dimension : 305mm(L) x 605mm(W) x 330mm(H)
- Net weight : 22KG



ACER BIOMEDICALS LLP

Office: Industrial Unit-80, Block-G, Sector-63, Noida, UP - 201301 (India)

Ph.: +91-120-4166921, E-mail: info@acerbiomed.com

Web.: www.acerbiomedicals.com, www.acerbiomed.com

Service Helpline
18008894232

