

# ABM

## ACER BIOMEDICALS

### Fingertip Pulse Oximeter FP-22



Accurate



Small and elegance



Reliable and stable



Weaving comfortable



Affinity design



#### Main Specification:

**Power Supply:** 2 AAA battery

**Display:** Color OLED

**Working current:** 30 mA

**Resolution:**  $\pm 2\%$

**Dimension:** 59\*36\*34mm

**Pulse rate(PR):** 30-250 bpm

**Blood oxygen saturation:** 35-100%

#### Main Features:

- Easy Operation
- Convenient detection with universal display
- Accurate blood oxygen saturation and pulse rate
- Anti-ambient light interference and anti-motion technology
- Low power consumption
- Battery display and low power indicator
- Auto shutdown if no signal

### ACER BIOMEDICALS LLP

Office: Industrial Unit-80, Block-G, Sector-63, Noida, UP - 201301 (India)  
Ph.: +91-120-4166921, E-mail: info@acerbiomed.com  
Web.: www.acerbiomedicals.com, www.acerbiomed.com

### Service Helpline

18008894232

