



PORTABLE SUCTION MACHINE 7E-H1



Technical Specification

Input Power	: 90 VA
Max vacuum	: ≥ 0.08 MPa
Max Air Flow	: 28 ± 2 L/min(without jar) 22 ± 2 L/min(with jar)
Operation Noise	: Approx, 60dB(A)(Without vacuum)
Liquid Collection Jar	: 1000ml, 1 piece
Operating Timew	: 30 mintes on, 30 MINTES off

ACER BIOMEDICALS LLP

Office: Industrial unit -80, Block-G, Sector-63, Noida, U.P.(India)
+91 120 4166921, E-mail: info@acerbiomed.com
Web: www.acerbiomedicals.com, www.acerbiomed.com
Toll free: 18008894232

Service Helpline
+91 120 4166921





PORTABLE SUCTION 7E-H1



TECHNICAL SPECIFICATIONS

Input power	:	90 VA
Max Vacuum	:	≥ 0.08 MPa
Max Air Flow	:	28 \pm 2 L/min (without jar) 22 \pm 2 L/min (with jar)
Operation noise	:	Approx. 60dB(A) (without vacuum)
Liquid Collection Jar	:	1000mL, 1piece
Operating time	:	30 minutes on, 30 minutes off

ACER BIOMEDICALS LLP

Corporate Office: D-10, Kamla Nagar, 1st Floor, New Delhi-110007,
Tel.: +91 11 4025 4357, 4025 4358,

E-mail: info@acerbiomed.com

Works: G-80, Sector-63, Noida, U.P.

Web: www.acerbiomedicals.com

Service Helpline
+ 91 20 4166921



