



ABJ200 Transcutaneous Jaundice Detector



Features:-

Advanced optics and software
simple to use
instant test results
Storage and memory function
No pain for infant

Specification:-

- 1) Display: LCD
- 2) Battery backup 1 horus (Recharable)
- 3) Power: AA 1.5Vx2 alkaline battery
- 4) Measurement range: 0.0mg/d L ~ 32.0 mg/d L
- 5) Measurement accuracy: Range is +1.5 mg/d L (25.5 pmol/L)
- 6) Charging preparation time: < 12 seconds
- 7) Record function: record 21 previous measurement values and review recorded data
- 8) Repeatability: < 10%

"Providing Clinical Excellence at the Point of Care"

ACER BIOMEDICALS LLP

Office: Industrial Unit-80, Block-G, Sector-63, Noida, UP - 201301 (India)
Ph.: +91-120-4166921, E-mail: info@acerbiomed.com
Web.: www.acerbiomedicals.com, www.acerbiomed.com

ServiceHelpline
18008894232

