



PORTABLE SUCTION MACHINE



7A-23D



Technical Specification

Input Voltage : AC220V, 50Hz

Max Negative Pressure: >80KPa

Operation Noise : < 60 dB(A)

Input Power : 180VA

Max Air Flow : 30±2 L/min(without jar)
26±2 L/min(with jar)

Liquid Collection Jar : 2500mL/pc, 2 Pieces

G.W./N.W : 5.2/13.5 Kgs

Size : 40.5x40x56 cm

ACER BIOMEDICALS LLP

Office: Industrial unit -80, Block-G, Sector-63, Noida, U.P.(India)

+91 120 4166921, E-mail: info@acerbiomed.com

Web: www.acerbiomedicals.com, www.acerbiomed.com

Toll free: 18008894232

Service Helpline
+91 120 4166921

