

## PORTABLE SUCTION MACHINE 7E-A



Oil-free  
pump



Over-flow  
protection

### Technical Specification

**Input Power** :90VA

**Max Vacuum** :0.75MPa

**Max Air Power** :  $\geq 17\text{L}/\text{min}$ (with jar)

**Operation Noise** : Approx.50db(A)(Without vacuum)(CE)  
Approx.50db(A)(Without vacuum)(Basic Sets)

**Liquid Collection Jar** :1000ml, 1piece

**Operating Time** : 30 minutes on, 30 minutes off

**G.W./N.W** :10/7.6 Kgs

**Size** : 54.5x37x31 cm

#### ACER BIOMEDICALS LLP

Office: Industrial unit -80, Block-G, Sector-63, Noida, U.P.(India)

+91 120 4166921, E-mail: [info@acerbiomed.com](mailto:info@acerbiomed.com)

Web: [www.acerbiomedicals.com](http://www.acerbiomedicals.com), [www.acerbiomed.com](http://www.acerbiomed.com)

Toll free: 18008894232

Service Helpline

+91 120 4166921

