




ECO Patient Monitor


Comprehensive and powerful monitoring solutions

Vital signs show perfectly

Suitable for ICU/CCU , Postoperative observatory, Emergency room, General ward, etc



Standard Interface



Max 8 waveform display simultaneously, max 6 leads ECG waveforms display

Product Feature


- Flexible, Light, Easy Operating
- Colorful and Clear 12.1" TFT screen, Backlight buttons available
- Various language operating systems
- Built-in detachable rechargeable lithium battery, 2.5hours battery life

Alarming – more safe and reliable


- 3 level Audible and visual alarming, Dual alarm light for physiological and technical alarming

Data Storage – stable and fast

- 400 groups NIBP measurements
- 72hours Trend data and Trend graphs
- 200 groups Alarm events reviewing
- 12hours Wave forms reviewing
- Supporting for Arrhythmia analysis and drug dose calculation



Big Font interface



Trend graphs

Standard Parameter: ECG, RESP, NIBP, SpO2, PR, TEMP
Optional: 2-IBP, EtCO2, Touch screen, Build in thermal printer



Durable probes connector with high quality



Build in Thermal Recorder (Optional)

ACER BIOMEDICALS LLP

Office: Industrial Unit-80, Block-G, Sector-63, Noida, UP - 201301 (India)
Ph.: +91-120-4166921, E-mail: info@acerbiomed.com
Web.: www.acerbiomedicals.com, www.acerbiomed.com

Service Helpline

18008894232

