

Pulse Oximeter P-Ultra



**8-inch
high-resolution
color TFT display**

**Lasting working
capability with
built-in rechargeable
lithium battery**

**Central monitoring
system applicable
to ICU/CCU/OT**

**Support solution:
wall mounting,
mobile trolley**

Features :

- 8 inch high denition color TFT display
- Standard parameter : SPO2 & PR
- 2-4 hours battery power supply
- Sound, light & text alarm, adjustable volume
- Connected to central monitoring System

SPO2:

Measurement range: 0-100%
Accuracy: 70-100%, $\pm 2\%$ 0-69%,
unspecified

Pulse rate:

Range 20-254bpm

Respiration:

Method RA-LL impedance
Rang

Weight: 1.700 kg

Electrical Specication:

100-240v AC, 50/60Hz

Max Input Power: 70VA

**Battery: 11.1V 4000mAh
rechargeable lithium battery**

**Standard configuration:
SpO2, PR**

ACER BIOMEDICALS LLP

Office: Industrial Unit-80, Block-G, Sector-63, Noida, UP - 201301 (India)
Ph.: +91-120-4166921, E-mail: info@acerbiomed.com
Web.: www.acerbiomedicals.com, www.acerbiomed.com

**Service Helpline
18008894232**

